

### **Amendments to the Claims:**

This listing of claims will replace all prior versions and listings of claims in the application:

### **Listing of Claims**

1-66. (Cancelled)

67. (Currently Amended) A quality control material that is not supplemented with bilirubin for use in a reagentless spectrophotometric apparatus, comprising, one or more substances that mimic ~~one or more analytes in whole blood, serum, plasma, urine, synovial fluid or cerebrospinal fluid~~ two or more analytes in a sample selected from the group consisting of plasma, serum, synovial fluid, cerebrospinal fluid, urine, mucus, lymphatic fluid, feces, semen, milk, cheese, cottage cheese, yogurt, ice cream, wine, beverages, semi-solid foods, soft solid foods, and any combination thereof ~~wherein said one or more analytes is selected from the group consisting of a simulator of turbidity, a perfluorocarbon like blood substitute, bilirubin, an indicator of hemolysis, a Hb-based blood substitute, methylene blue, met-Hb and biliverdin~~ , the quality control material being stabilized.

68. (Cancelled)

69. (Currently Amended) The quality control material of claim 68 ~~67~~, wherein ~~said one or more analyte~~ one of said two or more analytes is an indicator of hemolysis in ~~serum or plasma~~, wherein said indicator of hemolysis is one of total Hb, oxy-Hb or "total Hb minus met-Hb", ~~and wherein said quality control material is exposed to atmospheric conditions.~~

70. (Original) The quality control material of claim 69, wherein said one or more substances are selected from the group consisting of amaranth, acid fuchsin, basic

fuchsin, ponceau S, chromotrope 2R, phenol red, crystal ponceau, methyl orange, a Hb-based blood substitute, total Hb, oxy-Hb, carboxy-Hb, cyanmet-Hb, a polymer, and a protein.

71-72. (Cancelled)

73. (Currently Amended) A The quality control material of claim 67 wherein one of said two or more analytes is an indicator of hemolysis ~~for use in a reagentless spectrophotometric apparatus, comprising, one or more substances that mimics an~~ said indicator of hemolysis, ~~wherein said one or more substances are~~ is selected from the group consisting of total Hb, oxy-Hb, "total Hb minus met-Hb," cyanmet-Hb, amaranth, acid fuchsin, basic fuchsin, ponceau S, chromotrope 2R, phenol red, crystal ponceau, methyl orange, a Hb-based blood substitute, carboxy-Hb, a polymer, and a protein, ~~and wherein said quality control material is exposed to atmospheric conditions.~~

74-78. (Cancelled)

79. (Currently Amended) A The quality control material of claim 67 ~~for use in a reagentless spectrophotometric apparatus, comprising, one or more substances that mimics one or more of, wherein said two or more analytes area selected from the group consisting of an indicator of hemolysis,~~ biliverdin, bilirubin, methylene blue, met-Hb, a simulator of turbidity, ~~a perfluorocarbon like blood substitute~~ a blood substitute that is a milky-white emulsion and, a Hb-based blood substitute.

80. (Currently Amended) A The quality control material of claim 79 ~~for use in a reagentless spectrophotometric apparatus, comprising, one or more substances that mimics one or more of, an indicator of hemolysis, biliverdin, bilirubin, methelene blue, met Hb, a simulator of turbidity, a perfluorocarbon like blood substitute, a Hb-based blood substitute wherein said indicator of hemolysis is selected from the group consisting of total Hb, oxy-Hb and "total Hb minus met-Hb."~~

81. (Cancelled)

82. (Currently Amended) A The quality control material ~~for use in a reagentless spectrophotometric apparatus, comprising, of claim 67, wherein said~~ one or more substances ~~that~~ mimics an indicator of hemolysis, wherein said substance is selected from the group consisting of total-Hb and oxy-Hb, wherein said oxy-Hb accounts for about 95% of total Hb, or said total-Hb comprises about 95% oxy-Hb, and wherein said quality control material is exposed to atmospheric conditions.

83-97. (Cancelled)

98. (New) The quality control material of claim 67, wherein said one or more substances exhibits a negative slope in the absorbance spectrum for a continuous spectral segment of about 5nm to about 400nm, within the wavelength range of about 450nm to about 3000nm.

99. (New) The quality control material of claim 67, wherein said two or more analytes produce absorbances within the wavelength range of about 450nm to about 3000nm.

100. (New) The quality control material of claim 67, wherein said one or more substances are selected from the group consisting of methyl orange, phenol red, basic fuchsin, acid fuchsin, crystal ponceau, amaranth, ponceau S, chromotrope 2R, malachite green, brilliant green, toluidine blue O, evans blue, methylene blue, crystal violet, copper sulfate, total Hb, oxy-Hb, carboxy-Hb, cyanmet-Hb, a lipid emulsion, a blood substitute that is a milky-white emulsion, a Hb-based blood substitute, a polymer, a protein, and any combination thereof.

101. (New) The quality control material of claim 67, wherein the one or more substances do not mimic themselves.